

4. Quantification by GC/negative ion chemical ionization (NICI) MS is an extremely sensitive method to analyze F₂-IsoPs. Compounds are converted to pentafluorobenzyl (PFB) ester, trimethylsilyl (TMS) ether derivatives for this analysis. The lower limits of detection of F₂-IsoPs is in the range of 3–15 fmol using a deuterated standard with a blank of less than 5 parts per thousand. Thus, it is not necessary to assay more than 1–3 mL of a fluid such as plasma, 0.2 mL urine, or 50–1000 mg tissue.
5. Quantification of F₂-IsoPs in biological samples using the methods outlined requires a mass spectrometer with NICI capabilities. The use of NICI MS enhances the sensitivity of the assay by orders of magnitude compared to the use of electron ionization MS and electrospray ionization MS, although the use of capillary column technology will likely make electrospray ionization MS as sensitive in the future.

Acknowledgements

Supported by NIH grants DK48831, GM42056, GM15431, DK26657, CA77839, and CA68435. J.D.M. is the recipient of a Burroughs Wellcome Fund Physician Scientist Award in Translational Research.

References

1. Halliwell, B. and Gutteridge, J. M. C. (1990) Role of free radicals and catalytic metal ions in human disease: an overview. *Methods Enzymol.* **186**, 1–85.
2. DeZwart, L. L., Meerman, J. H. N., Commandeur, J. N. M., and Vermeulen, P. E. (1999) Biomarkers of free radical damage applications in experimental animals and in humans. *Free Rad. Biol. Med.* **26**, 202–226.
3. Halliwell, B. and Grootveld, M. (1987) The measurement of free radical reactions in humans. *FEBS Lett.* **213**, 9–14.
4. Morrow, J. D., Hill, K. E., Burk, R. F., Nammour, T. M., Badr, K. F., and Roberts, L. J. (1990) A series of prostaglandin F₂-like compounds are produced in vivo in humans by a non-cyclooxygenase, free radical catalyzed mechanism. *Proc. Natl. Acad. Sci. USA* **87**, 9383–9387.
5. Taber, D. F., Morrow, J. D., and Roberts, L. J. (1997) A nomenclature system for the isoprostanes. *Prostaglandins* **53**, 63–67.
6. Morrow, J. D. and Roberts, L. J. (1997) The isoprostanes: unique bioactive products of lipid peroxidation. *Prog. Lipid Res.* **36**, 1–21.
7. Morrow, J. D., Awad, J. A., Boss, H. J., Blair, I. A., and Roberts, L. J. (1992) Non-cyclooxygenase-derived prostanoids (F₂-isoprostanes) are formed in situ on phospholipids. *Proc. Natl. Acad. Sci. USA* **89**, 10,721–10,725.
8. Morrow, J. D. and Roberts, L. J. (1999) Mass spectrometric quantification of F₂-isoprostanes in biological fluids and tissues as measure of oxidant stress. *Methods Enzymol.* **300**, 3–12.

Formation of Apolipoprotein AI-AII Heterodimers by Oxidation of High-Density Lipoprotein

Audric S. Moses and Gordon A. Francis

1. Introduction

Uptake of oxidized low-density lipoprotein (LDL) by scavenger receptors allows cholesterol to be accumulated in an unregulated fashion by cells in the artery wall, an event considered pivotal to the formation of atherosclerotic plaques (1). High-density lipoprotein particles (HDL) are believed to protect against atherosclerosis in part by removing excess cholesterol from the artery wall (2). HDL, however, is itself susceptible to oxidation by numerous in vitro methods to an extent equal to, or greater than, LDL (reviewed in ref. 3). In vitro oxidation of HDL by peroxidase-generated tyrosyl radical ("tyrosylated HDL"), one mechanism of phagocyte-mediated oxidation (4), results in a markedly enhanced ability of HDL to promote the mobilization of cholesterol from a variety of cultured cells (5–7). If occurring in vivo, this oxidative modification of HDL may actually enhance its ability to protect against atherosclerosis.

Oxidation of HDL results in the covalent crosslinking of its two major apolipoproteins, apoAI and apoAII. Isolation of the monomeric and crosslinked apolipoprotein subspecies of tyrosylated HDL allowed us to identify the component responsible for enhanced cholesterol mobilization by these particles to be apolipoprotein AI-AII heterodimers (8). Separation of the individual tyrosylated HDL protein subspecies could not be accomplished using size exclusion or reverse-phase high-performance liquid chromatography (HPLC), or fast-performance liquid chromatography (FPLC). This separation was eventually achieved using preparative electrophoresis combined with elution of

From: *Methods in Molecular Biology*, vol. 188: *Oxidative Stress Biomarkers and Antioxidant Protocols*
 Edited by: D. Armstrong © Humana Press Inc., Totowa, NJ